

ASKIR C30 BR POLYCARBONATE JARS 2X2L SUCTION UNIT



Ref. BCASPI30



29/08/2024 14:49:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ASKIR C30 BR POLYCARBONATE JARS 2X2L SUCTION UNIT



Ref. BCASPI30

Intended for professional use.

Description

Ca-mi Askir C30 BR Suction Unit, autoclavable polycarbonate jars, 2 x 2000 ml

For professional use

Mobile surgical Suction Unit for draining body fluids and for oral, nasal and tracheal aspiration in adults and children.

Equipped with a powerful suction system that can drain up to 36 l/min max., two collection jars and a system on five casters, making it easy to transport.

Main applications: Emergency Department, Endoscopy, Gastroenterology, General Surgery, Obstetrics, Dentistry.

- Two different operating options: The AC/DC adapter and rechargeable batteries make this surgical Suction Unit an approved multipurpose unit, ensuring continuous interventions without the risk of overheating.

- Equipped with a visual and audible alarm system indicating if the battery level is low, the main system is also equipped with a pressure gauge (bar & kPa) and a vacuum regulator.

-

Accessories provided: 2 fluid collection jars with an anti-reflux valve with different options / 8 x 14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm / Antibacterial & hydrophobic filters (single patient) / CH20 cannula / Universal AC/DC adapter (100-240V 50/60Hz)-

Engine: Oil-free piston pump requiring no maintenance

- Power supply: 14V 4A with Universal AC / DC adapter, Rechargeable Internal Battery (12V 4A) / Autonomy 60min – Charging time 240min

-

Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg

-

Max free air flow: 36 l/min

-

Noise level: 65 dBA

-

Operating cycle: uninterrupted (connected to the mains)

-

Weight: 7.89 kg-

Dimensions: 32 x 30 x 99 cm

 29/08/2024 14:49:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ASKIR C30 BR POLYCARBONATE JARS 2X2L SUCTION UNIT



Ref. BCASPI30

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139166	HS Code	90189084
Warranty (months)	36	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	11.6
Bare Unit / Weight (kg)	7.89		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			

29/08/2024 14:49:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr