

ASKIR C30 BR ASKIR C30 BR SUCTION UNIT UU FLOVAC 2X2L



Ref. BCASPI32



29/08/2024 14:50:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ASKIR C30 BR ASKIR C30 BR SUCTION UNIT UU FLOVAC 2X2L



Ref. BCASPI32

Intended for professional use.

Description

Ca-mi Askir C30 BR Suction Unit, Flovac® single-use pouches, 2 x 2000 ml

Mobile surgical Suction Unit designed to drain body fluids and perform oral, nasal and tracheal aspiration in adults and children.

Equipped with two collection jars and a five-wheel system, the Askir C20 BR is easy to carry.

Main applications: Emergency Department, Endoscopy, Gastroenterology, General Surgery, Obstetrics, Dentistry.

- Two different operating options: The AC/DC adapter and rechargeable batteries make this surgical Suction Unit an approved multipurpose unit, ensuring continuous interventions without the risk of overheating.

- Equipped with a visual and audible alarm system indicating if the battery level is low, the main system is also equipped with a pressure gauge (bar & kPa) and a vacuum regulator.

-

2 Single Use Bags Flovac® 2000ml

-

Supplied accessories: Fluid collection jar with anti-reflux valve with different options / 8 x 14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm / Antibacterial & hydrophobic filter (single patient) / CH20 cannula / Power cord with Schuko plug and universal AC/DC adapter

- Engine: Oil-free piston pump requiring no maintenance

- Power supply: 14V 4A with Universal AC / DC adapter, Rechargeable Internal Battery (12V 4A) Autonomy 60min – Charging time 240min- Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg- Max free air flow: 36 l/min- Noise level: 65 dBA (with jar) / 68.5 dBA (without jar)- Operating cycle: uninterrupted (connected to the mains)

-

Weight: 7.89 kg

-

Dimensions: 32 x 30 x 99 cm

 29/08/2024 14:50:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ASKIR C30 BR ASKIR C30 BR SUCTION UNIT UU FLOVAC 2X2L



Ref. BCASPI32

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139180	HS Code	90189084
Warranty (months)	36	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	10
Bare Unit / Weight (kg)	7.89		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			

29/08/2024 14:50:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr