

## ASKIR C30 ASKIR C30 SUCTION UNIT UU FLOVAC 2X2L



Ref. BCASPI37



29/08/2024 14:52:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## ASKIR C30 ASKIR C30 SUCTION UNIT UU FLOVAC 2X2L



Ref. BCASPI37

Surgical Suction Unit on wheels and mobile, intended to drain body fluids and to practice oral, nasal and tracheal aspiration in adults and children.

### Description

Ca-mi Askir C30 Suction Unit, Flovac® single-use pouches, 2 x 2000 ml

Surgical Suction Unit on wheels and mobile, intended to drain body fluids and to practice oral, nasal and tracheal aspiration in adults and children.

Intended for professional use and equipped with a powerful suction system that can drain up to 40 l/min max, Askir C30 is equipped with two collection jars and a five-caster system making it easy to transport.

Available in Basic and FS

versions- BASIC FLOVAC version, including: Fluid collection jar with anti-reflux valve / 8 x 14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12 mm / Antibacterial & hydrophobic filter (single patient) / Power cord with Schuko plug

-

Flovac® single-use pockets-

Main

system equipped with a pressure gauge (bar & kPa) and a vacuum regulator

- Engine: Oil-free piston pump requiring no maintenance

-

Power supply: 230V-50Hz

-

Maximum vacuum (adjustable): -0.80 Bar -80 kPa -600 mmHg

-

Max free air flow: 40 l/min

-

Noise level: 61 dBA

- Weight: 6.2 kg / Dimensions: 32 x 99 x 30 cm

### Features



29/08/2024 14:52:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

|                             |     |                      |     |
|-----------------------------|-----|----------------------|-----|
| Packed unit / Length (cm)   | nc  | Packed unit / Weight | 9.4 |
| Bare Unit / Weight (kg)     | 6.5 |                      |     |
| SPECIFIC STORAGE CONDITIONS |     |                      |     |
| Minimum temperature         | -20 | Max temperature      | 45  |
| Minimum humidity            | NA  | Maximum humidity     | 85  |
| OTHER FEATURES              |     |                      |     |