

SUCTION UNIT CA-MI ASKIR C30 FS BOTTLE/ BAG UU FLOVAC 2X2L CDE/FOOT



Ref. BCASPI38



 29/08/2024 14:52:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCTION UNIT CA-MI ASKIR C30 FS BOTTLE/ BAG UU FLOVAC 2X2L CDE/FOOT



Ref. **BCASPI38**

Intended for the use of professionals.

Description

C30 Suction Unit, Flovac® single-use pouches, 2 x 2000 ml, Foot control-Foot

control

Surgical Suction Unit on wheels and mobile, intended to drain body fluids and to practice oral, nasal and tracheal aspiration in adults and children.

Intended for professional use and equipped with a powerful suction system that can drain up to 40 l/min max, Askir C30 is equipped with two collection jars and a system on five casters making it easy to transport.

- Main system equipped with a pressure gauge (bar & kPa) and a vacuum regulator
- Available in two versions: Basic or FS
- Accessories for the Basic version:

-

Fluid collection jar with an anti-reflux valve

-

8 x 14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12 mm

-

Antibacterial & hydrophobic filter (single patient)-

Power cord with Schuko plug

- Accessories for the FS version

:-

Fluid collection jar with an anti-reflux valve

-

8 x 14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12 mm

 29/08/2024 14:52:31 - For information purposes only, can be subject to change
Antibacterial & hydrophobic filter (single patient)

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

- technical characteristics



Engine: Oil-free piston pump requires no maintenance

3/3

Power supply: 230V-50Hz

Maximum vacuum (adjustable): -0.80 Bar -80 kPa -600 mmHg

SUCTION UNIT CA-MI ASKIR C30 FS BOTTLE/ BAG UU FLOVAC 2X2L CDE/FOOT



Max free air flow: 40 l/min

Noise level: 61 dBA

Ref. **BCASPI38**

Weight: 6.2 kg

Features

Dimensions: 32 x 99 x 30 cm

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139241	HS Code	90189084
Warranty (months)	36	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	9.15
Bare Unit / Weight (kg)	6.5		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			

📅 29/08/2024 14:52:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr