

SUCCION UNIT NEW HOSPIVAC 400, 2X2L POLYC JARS, VERSION FULL 2



Ref. BCASPI41



29/08/2024 14:55:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SUCCION UNIT NEW HOSPIVAC 400, 2X2L POLYC JARS, VERSION FULL 2



Ref. BCASPI41

Description

-with foot control and communicating jars

Designed to drain the patient's body fluids, gases, tissues and bones, during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

- The precise dashboard and the full range of accessories that come with it make it the ideal device for surgical aspiration and liposuction.
- Available in 3 versions: Basic, FS or full• Accessories for the Basic version:- Fluid collection jar with an anti-reflux valve

8 x 14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12 mm

- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko

plug• Accessories for the FS version:

- Fluid collection jar with an anti-reflux

valve- 8 x 14 Silicone pipes (autoclavable) and conical connectors ø 10-11-12 mm- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko plug

- Pedal

• Accessories for the Full version:

- Fluid collection jar with an anti-reflux valve

- 8 x 14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12 mm

- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko plug

- Pedal

- Electronic Transfer System between collection jars

• Technical characteristics

:- Engine: Oil-free piston pump requiring no maintenance

- Power supply: 230V-50Hz

- Maximum vacuum (adjustable): -0.90 Bar -90 kPa -675 mmHg

29/08/2024 14:55:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCCION UNIT NEW HOSPIVAC 400, 2X2L POLYC JARS, VERSION FULL 2



Ref. BCASPI41

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139272	HS Code	90189084
Warranty (months)	36	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	30	Packed unit / Height (cm)	92
Packed unit / Length (cm)	50	Packed unit / Weight	22
Bare Unit / Weight (kg)	20		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	1	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 14:55:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr