

SUCTION UNIT CA-MI NEW HOSPIVAC 400, FLOVAC® SINGLE-USE POUCHES, 2 X 2000 ML-VERSION FS 2



Ref. BCASPI64



 29/08/2024 15:05:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 400, FLOVAC® SINGLE-USE POUCHES, 2 X 2000 ML-VERSION FS 2



Ref. **BCASPI64**

Description

-with foot control

Designed to drain the patient's body fluids, gases, tissues and bones, during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

- The precise dashboard and the full range of accessories that come with it make it the ideal device for surgical aspiration and liposuction.
- Available in 3 versions: Basic, FS or full• Accessories for the Basic version:- Fluid collection jar with an anti-reflux valve

8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm

- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko

plug• Accessories for the FS version:

- Fluid collection jar with an anti-reflux

valve- 8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko plug

- Pedal

• Accessories for the Full version:

- Fluid collection jar with an anti-reflux valve

- 8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm

- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko plug,

- Pedal

- Electronic Transfer System between collection jars

• Technical characteristics

:- Engine: Oil-free piston pump requiring no maintenance

- Power supply: 230V-50Hz

29/08/2024 15:05:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 400, FLOVAC® SINGLE-USE POUCHES, 2 X 2000 ML-VERSION FS 2



Ref. [BCASPI64](#)

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139500	HS Code	90189084
Warranty (months)	36	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	0.8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	3
Packed unit / Weight	20	Bare Unit / Weight (kg)	20
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 15:05:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr