

## SUCTION UNIT CA-MI NEW HOSPIVAC 400 FLAC / SACH UU FLOVAC 2X3L FS 3



Ref. BCASPI67



29/08/2024 15:07:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## SUCTION UNIT CA-MI NEW HOSPIVAC 400 FLAC / SACH UU FLOVAC 2X3L FS 3



Ref. BCASPI67

### Description

FS 3 version with foot control

Designed to drain the patient's body fluids, gases, tissues and bones during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

- The precise dashboard and the full range of accessories that come with it make it the ideal device for surgical aspiration and liposuction.
- Available in 3 versions: Basic, FS or full• Accessories for the Basic version:- Fluid collection jar with an anti-reflux valve

8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm

Antibacterial & hydrophobic filter (single patient)

Power cord with Schuko

plug• Accessories for the FS version:

Fluid collection jar with an anti-reflux

valve- 8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm- Antibacterial & hydrophobic filter (single patient)

Power cord with Schuko plug

Pedal

• Accessories for the Full version:

Fluid collection jar with an anti-reflux valve

8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm

- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko plug,

Pedal

Electronic Transfer System between collection jars

• Technical characteristics

:- Engine: Oil-free piston pump requiring no maintenance

Power supply: 230V-50Hz

Maximum vacuum (adjustable): 0.90 Bar -90 kPa -675 mmHg

29/08/2024 15:07:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## SUCTION UNIT CA-MI NEW HOSPIVAC 400 FLAC / SACH UU FLOVAC 2X3L FS 3



Ref. BCASPI67

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139531	HS Code	90189084
Warranty (months)	36	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	20
Bare Unit / Weight (kg)	20		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 29/08/2024 15:07:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)