

## SUCTION UNIT CA-MI NEW HOSPIVAC 400 FLAC / SACH UU FLOVAC 2X3L FULL



Ref. BCASPI68

CE  
0123

DISPOSITIF MEDICAL  
CLASSE  
**IIa**  
MEDICAL DEVICE

GARANTIE  
**3**  
3 YEARS  
GUARANTEE



 29/08/2024 15:08:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## SUCTION UNIT CA-MI NEW HOSPIVAC 400 FLAC / SACH UU FLOVAC 2X3L FULL



Ref. **BCASPI68**

### Description

Full 3 version with foot control and communicating jars

Designed to drain the patient's body fluids, gases, tissues and bones, during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

- The precise dashboard and the full range of accessories that come with it make it the ideal device for surgical aspiration and liposuction.
- Available in 3 versions: Basic, FS or full• Accessories for the Basic version:- Fluid collection jar with an anti-reflux valve

8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm

- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko

plug• Accessories for the FS version:

- Fluid collection jar with an anti-reflux

valve- 8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko plug

- Pedal

• Accessories for the Full version:

- Fluid collection jar with an anti-reflux valve

- 8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm

- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko plug,

- Pedal

- Electronic Transfer System between collection jars

• Technical characteristics

:- Engine: Oil-free piston pump requiring no maintenance

- Power supply: 230V-50Hz

29/08/2024 15:08:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## SUCTION UNIT CA-MI NEW HOSPIVAC 400 FLAC / SACH UU FLOVAC 2X3L FULL



Ref. **BCASPI68**

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139548	HS Code	90189084
Warranty (months)	36	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	20
Bare Unit / Weight (kg)	20		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

29/08/2024 15:08:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)