

JARS POLI 2L WITH LID FOR ASPIRET, ASKIR 20,30



Ref. BCASPI80



29/08/2024 17:26:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

JARS POLI 2L WITH LID FOR ASPIRET, ASKIR 20,30



Ref. BCASPI80

Description

Makrolon 2L autoclavable Polycarbonate jar up to maximum 121°C

+ lid with anti-reflux valve

Features

CE Marking	CE0123	Medical device	No
Bar code	3661040169316	HS Code	90330090
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.2
Bare Unit / Weight (kg)	0.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:26:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr