

COVER FOR JAR 1000 ML AND 2000 ML

Ref. BCASPI84



 29/08/2024 17:27:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

COVER FOR JAR 1000 ML AND 2000 ML



Ref. BCASPI84

Features

CE Marking	NONE	Medical device	No
Bar code	3661040164786	HS Code	90330090
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.2
Bare Unit / Weight (kg)	0.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

 29/08/2024 17:27:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr