

ANTIBACTERIAL AND HYDROPHOBIC FILTER FOR ASPIRET ASKIR 20 / 30 EMIVAC




Ref. BCASPI96



 29/08/2024 15:36:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ANTIBACTERIAL AND HYDROPHOBIC FILTER FOR ASPIRET ASKIR 20 / 30 EMIVAC



Ref. **BCASPI96**

Description

- Anti-bacterial and hydrophobic filter
- Diam 64mm / Connector 8mm

Features

CE Marking	NONE	Medical device	No
Bar code	8054619101209	HS Code	90192000
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.05
Bare Unit / Weight (kg)	0.05		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 15:36:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr