

## CA-MI NEW ASPIRET SUCTION UNIT BOTTLE/BAG UU FLOVAC 2L



Ref. BCFLOV30



29/08/2024 21:21:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CA-MI NEW ASPIRET SUCTION UNIT BOTTLE/BAG UU FLOVAC 2L



Ref. BCFLOV30

Fixed electric surgical Suction Unit for draining body fluids and for performing oral, nasal and tracheal aspiration in adults and children.

### Description

New Aspiret CA-MI Phlegm Vacuum - FLOVAC® 2L Single Use Bottle/Bag

Fixed electric surgical Suction Unit designed to drain body fluids and perform oral, nasal and tracheal aspiration in adults and children.

Suitable for tracheostomy patients and is ideal for minor surgical applications as well as postoperative home treatments.

Main applications: ENT, Home care, Bedroom, Tracheostomy

Composed of: New Aspiret Suction Unit + Flovac 2L Kit (1x Reusable Flovac 2L Jar / 1x Sterile UU Pocket 2L Flovac/ 1x Set tubes Ø6X10mm + connectors / 1x Round Discer for Flovac System / 1x Yellow Male Connector FLOVAC System (5 pcs) / 1x Configuration Manual Suction Units )

- Optional

roller foot- The main system is equipped with a pressure gauge (bar & kPa) and a vacuum regulator

-

Supplied accessories: Fluid collection jar with an anti-reflux valve with different options / 6 x 10 Silicone hoses (autoclavable) and conical connectors ø 8-9-10 mm / Antibacterial & hydrophobic filter (single patient) / Cannula CH20 / Power cord with Schuko plug- Motor:

Maintenance-free oil-free piston pump

-

Power supply: 230V-50Hz - other configurations available on request

-

Maximum vacuum (adjustable): 0.75 Bar -75 kPa -563 mmHg

-

Max free air flow: 15 l/min

-

Noise level: 59.6 dBA

- Operating

e: Uninterrupted  
29/05/2024 21:21:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CA-MI NEW ASPIRET SUCTION UNIT BOTTLE/BAG UU FLOVAC 2L



Ref. BCFLOV30

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040197227	HS Code	90189084
Warranty (months)	24	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	2.2
Bare Unit / Weight (kg)	2.2		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			

 29/08/2024 21:21:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)