

CA-MI NEW ASKIR 30 PROXIMITY ELECTR SUCTION UNIT, 2L JAR IN MAKROLON



Ref. BCPOLY401



📅 29/08/2024 21:17:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

CA-MI NEW ASKIR 30 PROXIMITY ELECTR SUCTION UNIT, 2L JAR IN MAKROLON



Ref. BCPOLY401

Fixed electric surgical Suction Unit, intended to drain body fluids and to perform oral, nasal and tracheal aspiration in adults and children.

Description

Suitable for tracheostomy patients, ideal for minor surgical applications as well as postoperative home treatments. The "PROXIMITY" function is used to activate or deactivate the device without touch contacts and prevents cross-contamination between two patients.

Main applications: Endoscopy, Gastroenterology, General surgery, Home care, Obstetrics, Room, Veterinary Med.

- Oil-free piston pump requiring no maintenance

- Main system is equipped with a pressure gauge (bar & kPa) and a vacuum regulator

- Supplied accessories: Fluid collection jar with anti-reflux valve with different options / 6 x 10 Silicone hoses (autoclavable) and tapered connectors ø 8-9-10 mm / Antibacterial & hydrophobic filter (single patient) / CH20 cannula / Power cord with Schuko

plug- Power supply: 230V-50Hz - other configurations available on request

- Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg

- Max free air flow: 40 l/min

- Noise level: 60.5 dBA

- Uninterrupted operating cycle

- Autoclavable Makrolon polycarbonate jar

- Flow regulation button and vacuum gauge

- Easily washable ABS shell


- Antibacterial filter

- Dimensions: 35 x 21 x 18 cm / Weight: 3kg60

 29/08/2024 21:17:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

CA-MI NEW ASKIR 30 PROXIMITY ELECTR SUCTION UNIT, 2L JAR IN MAKROLON



Ref. [BCPOLY401](#)

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040200149	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	35
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	18
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	21
Packed unit / Weight	3.6	Bare Unit / Weight (kg)	3.6
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 21:17:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr