

FOLDING SIDEBAR AND DISPLAY HOLDER




Ref. EKMWA04

 29/08/2024 17:20:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

FOLDING SIDEBAR AND DISPLAY HOLDER



Ref. **EKMWA04**

Description

Ergonomic restraint bar KERN MWA-A01 (with a foot for ergonomic mounting of the analytical device), which offers the patient maximum safety to climb on the weighing platform

For any measuring device requiring medical Class III approval, the accessories must be ordered at the same time as it so that the approval can be valid.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040173672	HS Code	84231090
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	40
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	101.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	56
Packed unit / Weight	nc	Bare Unit / Weight (kg)	4
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:20:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr