

MIXTER-O'SHAUGNESSY ARTERY FORCEPS 90° 16 CM



Ref. IP49016



29/08/2024 15:18:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MIXTER-O'SHAUGNESSY ARTERY FORCEPS 90° 16 CM



Ref. IP49016

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040138305	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	1
Packed unit / Length (cm)	22	Packed unit / Weight	0.07
Bare Unit / Weight (kg)	0.07		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 15:18:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr