

TROLLEY WITH WHEELS ELITE BAGS




Ref. MECHAR02



29/08/2024 17:32:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

TROLLEY WITH WHEELS ELITE BAGS



Ref. **MECHAR02**

External trolley, foldable, robust and lightweight.

Description

It folds quickly thanks to a button located on the handle. It has an elastic rope that allows you to completely fix the load.

Materials: steel tubes, ABS handle, PP base and PVC wheels

Wheel dimensions: Ø 77 x 22 mm

Color: black

Weight: 1.35 kg- Outdoor measures: 23 x

96.5 x 32 cm (folded: 23 x 44 cm)

Recommended max. capacity: 30 kg

- Can be used with models: Call, Podia, Doctor, Sport, Medic, Cool, Practi, Vet, Urban and Doctor

Features

CE Marking	NONE	Medical device	No
Bar code	8436564540621	HS Code	732690
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	25	Bare Unit / Width (cm)	23
Packed unit / Height (cm)	11	Bare Unit / Height (cm)	32
Packed unit / Length (cm)	48	Bare Unit / Length (cm)	96.5
Packed unit / Weight	1.8	Bare Unit / Weight (kg)	1.35
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:32:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr