

PLASTIC TROLLEY, 2 TRAYS 48,5 X 36,5 CM

Ref. OM483602




CE



 29/08/2024 21:15:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

PLASTIC TROLLEY, 2 TRAYS 48,5 X 36,5 CM

Ref. **OM483602**

Lightweight, handy and very robust pedestal table

Description

The quality of the casters and ergonomic handles ensure maximum comfort when travelling. - External dimensions: 48.5 x 36.5 x

- 77 cm, including handles
- Dimensions of the trays: 48.5 x 36.5 cm
- Maximum load: 100 kg per trolley and 60 kg per level
- Color: Light gray
- Materials: Aluminum upright, PP top
- 4 PVC swivel casters, Ø 75 mm including 2 brake
- Weight: 4 kg

Features

CE Marking	CE	Medical device	No
Bar code	3760206421302	HS Code	8716800000
Warranty (months)	12		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	41	Bare Unit / Width (cm)	36.5
Packed unit / Height (cm)	19	Bare Unit / Height (cm)	77
Packed unit / Length (cm)	54	Bare Unit / Length (cm)	48.5
Packed unit / Weight	7	Bare Unit / Weight (kg)	4
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	35
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 29/08/2024 21:15:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr