

ELECTRODES FOR TRAINING DEFIBRILLATOR (5 PAIRS)

Saver One®

Ref. SJACCE69



 29/08/2024 14:30:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ELECTRODES FOR TRAINING DEFIBRILLATOR (5 PAIRS)



Ref. **SJACCE69**

Description

Designed solely for training in the use of a defibrillator in cardiac resuscitation. Do not use on a human.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040137285	HS Code	90181990
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	16
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	22
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	23
Packed unit / Weight	0.1	Bare Unit / Weight (kg)	0.1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 14:30:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr