

## DISPOSABLE FLOVAC LINER 2L FOR ASPIRET, ASKIR 20, 30, C30BR



Ref. BCASPI101



29/08/2024 16:09:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DISPOSABLE FLOVAC LINER 2L FOR ASPIRET, ASKIR 20, 30, C30BR



Ref. BCASPI101

### Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	0000000000604	HS Code	90189084
Warranty (months)	12	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.2
Bare Unit / Weight (kg)	0.2		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

29/08/2024 16:09:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)