

CONTAINER FLOVAC 2L (REUSABLE) FOR ASPIRET, ASKIR 20, 30 AND C30BR



Ref. BCASPI104



 29/08/2024 16:09:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

CONTAINER FLOVAC 2L (REUSABLE) FOR ASPIRET, ASKIR 20, 30 AND C30BR



Ref. BCASPI104

Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	0000000000611	HS Code	4202121900
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	14
Packed unit / Length (cm)	34	Packed unit / Weight	0.366
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 16:09:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr