

SET OF SILICONE TUBES Ø8X14 JAR + FILTER FOR ASKIR C30BR




Ref. BCASPI116



29/08/2024 16:04:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SET OF SILICONE TUBES Ø8X14 JAR + FILTER FOR ASKIR C30BR



Ref. BCASPI116

Description

Set of pipes containing:- 1 pipe x diam 8 / ext 14mm x length

150

cm-1 pipe x diam 8 / ext 14mm x length 17 cm

Conical connectors 10-11-12 mm

Features

CE Marking	CE0123	Medical device	No
Bar code	3661040166346	HS Code	90330090
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.242
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 16:04:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr