

JARS IN POLICARBONATE (MAX 121°C) 400 ML WITH COVER FOR EMIVAC



Ref. [BCASPI77](#)

 29/08/2024 15:43:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

JARS IN POLICARBONATE (MAX 121°C) 400 ML WITH COVER FOR EMIVAC



Ref. [BCASPI77](#)

Features

CE Marking	NONE	Medical device	No
Bar code	3661040151076	HS Code	4202121900
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	18.4
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	9.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	27
Packed unit / Weight	0.15	Bare Unit / Weight (kg)	0.1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 15:43:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr