

JARS POLI 1L WITH LID FOR ASPIRET, ASKIR 20,30




Ref. BCASPI78



 29/08/2024 15:34:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

JARS POLI 1L WITH LID FOR ASPIRET, ASKIR 20,30



Ref. BCASPI78

Description

1L jar in Makrolon 2858 autoclavable at maximum 121 ° C

+ lid with an anti-reflux valve

Features

CE Marking	CE0124	Medical device	No
Bar code	3661040151083	HS Code	90330090
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	12
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	12
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	21
Packed unit / Weight	0.2	Bare Unit / Weight (kg)	0.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 15:34:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr