

SPECIAL DRESS HANDKERCHIEF ANALYZES ROW ELITE BAGS, BLUE



Ref. MEROW02_B



29/08/2024 16:44:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SPECIAL DRESS HANDKERCHIEF ANALYZES ROW ELITE BAGS, BLUE



Ref. MEROW02_B

Small insulated case very functional and handy, designed for home visits

Description

- It allows the professional to transport in good order and in optimal conditions the samples obtained- 1 small isothermal spongy support for various tube diameters: 3 large, 28 medium and 5 small- 1 hollow holder for syringes, elastics, needles and gloves
- 1 net bag with zipper for the cold accumulator and for common utensils- 2 external pockets for the container of bio-contaminated material pocket and for the glucometer- 1 hermetic pouch
- 1 absorbent sponge base
- Material: 340D fabric
- Weight: 480 g
- Outdoor measures: 26 x 12 x 12 cm / 3.7 L
- Maximum recommended capacity: 2 kg
- Certified for the transport of biological substances UN 3373.
- Sold empty

Features

CE Marking	NONE	Medical device	No
Bar code	8436564540348	HS Code	42029270
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	29	Bare Unit / Width (cm)	12
Packed unit / Height (cm)	11	Bare Unit / Height (cm)	12
Packed unit / Length (cm)	35	Bare Unit / Length (cm)	26
Packed unit / Weight	0.44	Bare Unit / Weight (kg)	0.48
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr