

## OSCILLATING PLASTER ELECTRIC SAW - VARIATOR - (MODEL GOLD +)

Ref. ZSE405602



27/08/2024 07:28:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## OSCILLATING PLASTER ELECTRIC SAW - VARIATOR - (MODEL GOLD +)



Ref. ZSE405602

### Description

Delivered with:- 1 saw blade 50 mm resin- 1 saw blade

65 mm resin- 2 keys

- Engine cooling by air circulation system
- Rotation speed control for optimal level of oscillation of saw blades
- Ergonomic switch
- Easy blade replacement- Electrical connection 220-240 Volt, 50/60 Hz.
- Rated power 250 Watt
- Oscillations ca. 12000 - 21000 min
- Weight 1.2 kg
- Protection category II
- Protection (IP) IP2X
- Sound level 75 dB. The noise level when working can exceed 85 dB. From 85 dB, it is advisable to wear hearing protection.
- Hand-arm vibration Less than 2.5 m/s.

### Features

CE Marking	CE	Medical device	Class I
Bar code	8600131047344	HS Code	84672290
Warranty (months)	24	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	20	Bare Unit / Width (cm)	20
Packed unit / Height (cm)	9	Bare Unit / Height (cm)	9
Packed unit / Length (cm)	44	Bare Unit / Length (cm)	44.5
Packed unit / Weight	2.38	Bare Unit / Weight (kg)	2.1
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	5	Max temperature	40
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

27/08/2024 07:28:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)