

## LID FOR BEDPAN (STAINLESS STEEL)

Ref. AI14833



29/08/2024 12:40:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LID FOR BEDPAN (STAINLESS STEEL)



Ref. AI14833

### Features

|                                    |               |                           |          |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking                         | CE            | Medical device            | Class I  |
| Bar code                           | 3661040117379 | HS Code                   | 90189090 |
| Warranty (months)                  | 12            | Country                   | Pakistan |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |          |
| Packed unit / Width (cm)           | 31            | Packed unit / Height (cm) | 2        |
| Packed unit / Length (cm)          | 31            | Packed unit / Weight      | 0.416    |
| Bare Unit / Weight (kg)            | 0.416         |                           |          |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |          |
| Minimum temperature                | NA            | Max temperature           | NA       |
| Minimum humidity                   | NA            | Maximum humidity          | NA       |
| <b>OTHER FEATURES</b>              |               |                           |          |
| Packed unit / Quantity             | 1             |                           |          |

29/08/2024 12:40:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)