

INTERFACE CABLE RS 232 FOR KERN MWS, MXS, MTS, MPS

KERN[®]


Ref. EKMP508



 27/08/2024 08:24:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

INTERFACE CABLE RS 232 FOR KERN MWS, MXS, MTS, MPS



Ref. EKMP508

Description

Kern RS 232C interface cable for connection to an external device

Features

CE Marking	NONE	Medical device	No
Bar code	4045761096950	HS Code	84231090
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 08:24:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr