

FILTER FOR NEW HOSPIVAC 400



Ref. BCASPI106

 29/08/2024 16:59:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

FILTER FOR NEW HOSPIVAC 400



Ref. BCASPI106

Description

Filter for Suction Unit HOSPIVAC 400, diam 90 with connector 11mm

Features

CE Marking	NONE	Medical device	No
Bar code	8054619103708	HS Code	90192000
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	6
Packed unit / Length (cm)	9	Packed unit / Weight	0.042
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 16:59:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr