

## SILICONE TUBES 8X14MMX150CM/8X14MMX17+CONNECTORS FOR JARS HOSPIVAC




Ref. BCASPI139



 29/08/2024 17:46:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## SILICONE TUBES 8X14MMX150CM/8X14MMX17+CONNECTORS FOR JARS HOSPIVAC



Ref. **BCASPI139**

### Description

Set of pipes containing:- 1 pipe x diam 8 / ext 14mm x length

150

cm-1 pipe x diam 8 / ext 14mm x length 17 cm

-  
Conical connectors 10-11-12 mm

### Features

CE Marking	CE0123	Medical device	No
Bar code	0000000000857	HS Code	90330090
Warranty (months)	12	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	18	Packed unit / Height (cm)	3
Packed unit / Length (cm)	17	Packed unit / Weight	0.25
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

29/08/2024 17:46:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)