

LEAD BATTERY FOR ASKIR 230/12 BR, 36 BR, C30 BR



Ref. BCASPI165



 29/08/2024 17:44:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LEAD BATTERY FOR ASKIR 230/12 BR, 36 BR, C30 BR



Ref. BCASPI165

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	3661040169095	HS Code	85065090
Warranty (months)	36	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:44:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr