

DISPOSABLE PROBE COVERS FOR THERMOSCAN PRO 6000, BOX OF 800

WelchAllyn®

Advancing Frontline Care™

Ref. DWPRO6K02



29/08/2024 17:57:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DISPOSABLE PROBE COVERS FOR THERMOSCAN PRO 6000, BOX OF 800

WelchAllyn®

Advancing Frontline Care™

Ref. **DWPRO6K02**

Description

This covers probe helps protect your ThermoScan Pro 6000 ear thermometer.

Single use for better efficiency and impeccable hygiene, this protective cap is perfectly suited to the Braun thermometer.

Its shape is studied to suit children and babies.

Colour: transparent

Features

CE Marking	CE	Medical device	Class I
Bar code	00732094321043	HS Code	9025900090
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	8
Packed unit / Length (cm)	30	Packed unit / Weight	0.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 29/08/2024 17:57:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr