

BATTERY BOX 2 PIN CONNECTOR FOR MCC, MCB

KERN[®]

Ref. EKACCE22



 29/08/2024 18:16:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

BATTERY BOX 2 PIN CONNECTOR FOR MCC, MCB



Ref. **EKACCE22**

Description

2-pin connector battery case for MCC

Features

CE Marking	NONE	Medical device	No
Bar code	3661040173481	HS Code	84231090
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 18:16:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr