

ELITE BAGS - PRACTI - BLUE



Ref. MEPRACT02_N



29/08/2024 19:33:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ELITE BAGS - PRACTI - BLUE



Ref. MEPRACT02_N

Health assistance bag for home consultations or emergency.

Description

Its excellent internal distribution makes it possible to immediately identify materials.

- 2 extractable compartments with transparent front part

- Isothermal amp maker

- External dividers

- Inner and outer pockets- Measures: 40 x

30 x 16 cm

- Capacity: 19.20 L

- Weight: 1.82 kg

- Material: 420D polyester

- Washable

- Maximum recommended capacity: 4 kg

- Sold empty

Features

CE Marking	NONE	Medical device	No
Bar code	8436564541819	HS Code	42029900
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	14	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	31	Bare Unit / Height (cm)	16
Packed unit / Length (cm)	43	Bare Unit / Length (cm)	40
Packed unit / Weight	2	Bare Unit / Weight (kg)	1.82

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Black
--------	-------

29/08/2024 19:33:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

+33(0)4 42 90 31 31

info@gsh-med.fr

www.gsh-med.fr