

FOOT STEP - 2 STEPS - EPOXY

Ref. OS0052E



29/08/2024 18:33:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

FOOT STEP - 2 STEPS - EPOXY



Ref. OS0052E

Description

- Structure in epoxy lacquered steel tube
- Steps covered with non-slip plastic
- Maximum load: 150 kg- Dimensions: (L)40 x (W)24.5 x (H)34 cm
- Height of the first step: 16.5 cm

Stainless steel or Epoxy: what difference?

Epoxy: the economic choice. The steel is simply coated with a resin to delay corrosion.

Stainless steel: it is steel modified to become completely STAINLESS steel, it does not rust.

And did you know?

Today, more than 90% of stainless steel is made from recycled steel, it's good for the environment!

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040168104	HS Code	94029000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	35	Bare Unit / Width (cm)	24.5
Packed unit / Height (cm)	12	Bare Unit / Height (cm)	34
Packed unit / Length (cm)	57	Bare Unit / Length (cm)	40
Packed unit / Weight	6	Bare Unit / Weight (kg)	5.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 18:33:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr