

STAINLESS STEEL SCREEN 4 PANELS WITH STRAINED WHITE CURTAINS- HOLTEX



Ref. OS0064



29/08/2024 18:34:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

STAINLESS STEEL SCREEN 4 PANELS WITH STRAINED WHITE CURTAINS- HOLTEX



Ref. OS0064

Description

- Panels (155 x 61 cm) with non-separable casters without brakes
- Chassis made of stainless steel tube
- Washable plastic curtains
- Fire resistance classification: M3
- Double casters. Diameter: 60 mm
- Height: 175 cm

Features

CE Marking	NONE	Medical device	No
Bar code	3661040168111	HS Code	9403899000
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	64	Packed unit / Height (cm)	10
Packed unit / Length (cm)	177	Packed unit / Weight	13
Bare Unit / Weight (kg)	10.6		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 18:34:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr