

PEDAL BIN WITH LID, 20 LITERS S/S HOLTEX

Ref. OSA0101



29/08/2024 18:24:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PEDAL BIN WITH LID, 20 LITERS S/S HOLTEx



Ref. OSA0101

Description

- Stainless steel- Closed height: 43 cm
- Open lid height: 74 cm
- Outer diameter: 29 cm
- Removable plastic inner bucket

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040168241	HS Code	7323990000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	32	Bare Unit / Width (cm)	29
Packed unit / Height (cm)	45	Bare Unit / Height (cm)	43
Packed unit / Length (cm)	32	Packed unit / Weight	3
Bare Unit / Weight (kg)	2.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 18:24:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr