

## INSTRUMENT TROLLEY S.S 3 SHELVES WITH SIDE RAILS 60X40X80 CM



Ref. OSA6041



29/08/2024 18:26:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## INSTRUMENT TROLLEY S.S 3 SHELVES WITH SIDE RAILS 60X40X80 CM



Ref. OSA6041

High quality pedestal table, perfect for storing the necessary surgical instruments or consumables

### Description

- Round tube structure in non-magnetic stainless steel for a better grip
- Stainless steel trays
- Removable galleries- Silent rubber wheels
- Permissible load: 80 kg
- Wheel diameter: 75 mm with brakes

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040168296	HS Code	9402900000
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	45	Bare Unit / Width (cm)	40
Packed unit / Height (cm)	14	Bare Unit / Height (cm)	80
Packed unit / Length (cm)	85	Bare Unit / Length (cm)	60
Packed unit / Weight	12.2	Bare Unit / Weight (kg)	12.2
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

29/08/2024 18:26:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)