


HEAD IMMOBILIZER, ORANGE

Ref. SBRAN10_O



29/08/2024 16:52:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

HEAD IMMOBILIZER, ORANGE



Ref. **SBRAN10_O**

Description

- Manufactured with high density plastic materials, wide holes at the ears to monitor the auditory canal of the patient
- It can be used with a stretcher and a hard board
- Minimum Interference with radiography systems or MRI.
- Plastic Coated, water-resistant and easy to clean
- Dimensions (L x W x H): 40 x 25 x 16 cm
- Weight: 1 kg
- Single colour: orange

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040094656	HS Code	94029000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	26	Bare Unit / Width (cm)	25
Packed unit / Height (cm)	14	Bare Unit / Height (cm)	16
Packed unit / Length (cm)	43	Bare Unit / Length (cm)	40
Packed unit / Weight	0.95	Bare Unit / Weight (kg)	1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 16:52:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr