

## PE/PLASTIC SCOOP STRETCHER, YELLOW

Ref. SBRAN13\_J



29/08/2024 16:53:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## PE/PLASTIC SCOOP STRETCHER, YELLOW



Ref. **SBRAN13\_J**

=====

### Description

- «Spoon» function performed by the two longitudinal blades, to position the stretcher under the patient without moving him. The scoop stretcher can be used to lift, move, or transport patients.
- Length adjustable on 4 levels according to the patient's size, from 162 to 234 cm
- High resistance aluminium structure
- Dimensions (L x W x H): 162-234 x 44 x 7 cm
- Weight: 7 kg
- Maximum load carried: 220 kg

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040094670	HS Code	94029000
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	47	Bare Unit / Width (cm)	45
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	7
Packed unit / Length (cm)	165	Bare Unit / Length (cm)	162
Packed unit / Weight	10.8	Bare Unit / Weight (kg)	8
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

29/08/2024 16:53:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)