

ALUMINIUM ALLOY STAIR STRETCHER - W 38 CM - 2 WHEELS, ORANGE



Ref. SBRAN14_O



29/08/2024 16:54:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ALUMINIUM ALLOY STAIR STRETCHER - W 38 CM - 2 WHEELS, ORANGE



Ref. **SBRAN14_O**

Ideal for transferring a patient, going up and down stairs

Description

- High strength aluminium alloy structure
- Ideal for transferring a patient, climbing up and down the stairs
- Two fixed rear wheels allow movement on most surfaces with minimal effort
- The safety of patients during the transport process is ensured by 3 restraining straps with quick release buckle
- At the front of the stretcher, two retractable arms
- Dimensions (L x W x H): 62 X 51 x 91 cm, unfolded; 90 x 17 x 59 cm, folded
- Weight: 8 kg
- Maximum load carried: 159 kg
- Soft seat
- Seat height: 49 cm
- Single colour: orange

Features

| | | | |
|------------------------------------|---------------|-------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040094687 | HS Code | 94029000 |
| Warranty (months) | 12 | Country | China |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 56 | Bare Unit / Width (cm) | 52 |
| Packed unit / Height (cm) | 16 | Bare Unit / Height (cm) | 92 |
| Packed unit / Length (cm) | 98 | Bare Unit / Length (cm) | 84 |
| Packed unit / Weight | 10 | Bare Unit / Weight (kg) | 9 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

+33(0)4 42 90 31 31

info@gsh-med.fr

www.gsh-med.fr