

PROBE FOR VASCULAR DOPPLER / 5 MHZ



Ref. 743311



 30/08/2024 21:32:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

PROBE FOR VASCULAR DOPPLER / 5 MHZ



Ref. 743311

Description

5 MHz probe - vascular probe for doppler

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446008178	HS Code	90181990
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.145
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	55
Minimum humidity	5	Maximum humidity	93
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 21:32:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr