

## SUNTECH CT50 BP STANDARD SPO2

# SUNTECH

Ref. DSU222201



CE  
0413

DISPOSITIF MEDICAL  
CLASSE  
**IIb**  
MEDICAL DEVICE

GARANTIE  
2 ANS  
**2**  
2 YEARS  
GUARANTEE



30/08/2024 09:27:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

# SUNTECH CT50 BP STANDARD SPO2



Ref. **DSU222201**

**The world reference for the measurement of vital parameters. The CT 50 is a compendium of state-of-the-art technology**

## Description

The SUNTECH CT 50 monitor is a condensed version of the latest technology.  
 The world reference for measuring vital signs.  
 Nellcor SpO2 technology, optional HIS, USPx4, Ethernet LAN and WIFI connectivity.  
 3 blood pressure measurement modes: Automatic, Auscultatory and average mode.  
 8-inch colour touch screen, optional barcode reader, optional printer and optional wheeled base

## Features

CE Marking	CE0413	Medical device	Class IIb
Bar code	3700446025113	HS Code	9018199560
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	4.9
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 30/08/2024 09:27:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
 13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)