

SUNTECH CT50 BP STANDARD SPO2 TEMP PRINTER

SUNTECH

Ref. DSU222204

CE
0413

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
2
ANS
2 YEARS
GUARANTEE



📅 30/08/2024 10:52:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

SUNTECH CT50 BP STANDARD SPO2 TEMP PRINTER



Ref. **DSU222204**

The world reference for the measurement of vital parameters. The CT 50 is a compendium of state-of-the-art technology

Description

The SUNTECH CT 50 monitor is a compendium of state-of-the-art technology. The world reference for the measurement of vital signs. Nellcor SpO2 technology, optional HIS, USPx4, Ethernet LAN and WIFI connectivity. 3 blood pressure measurement modes: Automatic, Auscultatory and medium mode. 8-inch colour touch screen, optional barcode reader, optional printer and optional wheel stand

Features

CE Marking	CE0413	Medical device	Class IIb
Bar code	3661040214238	HS Code	9018199560
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	32	Packed unit / Height (cm)	33
Packed unit / Length (cm)	37	Packed unit / Weight	5.644
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 10:52:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr