

CABLE TANGO RS232 INTERFACE FOR CASE 8000 BOX 16...

SUNTECH

Ref. DSU230748

CE
0413



 30/08/2024 10:15:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

CABLE TANGO RS232 INTERFACE FOR CASE 8000 BOX 16...

SUNTECH

Ref. DSU230748

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	3700446012410	HS Code	9018110000
Warranty (months)	3	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	20	Packed unit / Height (cm)	5
Packed unit / Length (cm)	22	Packed unit / Weight	0.12
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 10:15:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr