

USB PC LINK CABLE - OSCAR 2 MAPA 24 H

SUNTECH

Ref. DSU230912



CE
0413



 30/08/2024 10:06:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

USB PC LINK CABLE - OSCAR 2 MAPA 24 H

SUNTECH

Ref. DSU230912

Description

Oscar 2TM cables.
Cables to connect Oscar 2TM to a computer.
To program and import patients' data.


Features

CE Marking	CE0413	Medical device	Class I
HS Code	9018199560	Warranty (months)	3
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.07
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 10:06:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr