

SET OF 5 SUNTECH MONOBLOC DOUBLE TUBE REUSABLE ARMBANDS SMALL ADULT (17-25CM) - SCREW CONNECTOR

SUNTECH

Ref. DSU231034

CE



 30/08/2024 09:47:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SET OF 5 SUNTECH MONOBLOC DOUBLE TUBE REUSABLE ARMBANDS SMALL ADULT (17-25CM) - SCREW CONNECTOR



Ref. [DSU231034](#)

Description

Set of 5 Suntech turquoise nylon armbands.
Double tubing and screw connection, small adult size 17-25cm

Features

CE Marking	CE	Medical device	Class I
Bar code	0120840935104084	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.348
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 09:47:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr