

## TIGHT 38-50CM X20

Ref. DSU231069

# SUNTECH

CE



 30/08/2024 09:51:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## TIGHT 38-50CM X20

SUNTECH

Ref. **DSU231069**

### Description

Set of 20 Suntech "SOFT" single-patient armbands, brown.  
Bib shorts size 38-50cm, rectus bayonnette connection.

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014407	HS Code	9018905040
Country	China		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	2.49
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 30/08/2024 09:51:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)