


## SPO2 REUSABLE DIGITAL SENSOR FOR ADULT 3 METERS

SUNTECH

Ref. DSU250899

CE  
0413


DISPOSITIF MEDICAL  
CLASSE  
**II<sub>b</sub>**  
MEDICAL DEVICE



 30/08/2024 09:26:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

# SPO2 REUSABLE DIGITAL SENSOR FOR ADULT 3 METERS

SUNTECH

Ref. DSU250899

## Description

SpO2 articulated finger sensor for adult cable with a length of 3 meters for Suntech CT 40 station


## Features

CE Marking	CE0413	Medical device	Class IIb
Bar code	0108698870322340	HS Code	90181990
Warranty (months)	3	Country	Turkey
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	3
Packed unit / Length (cm)	13	Packed unit / Weight	0.115
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 30/08/2024 09:26:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)