

SPO2 REUSABLE DIGITAL SENSOR FOR ADULT

SUNTECH

Ref. DSU250900

CE
0413

DISPOSITIF MEDICAL
CLASSE
II_b
MEDICAL DEVICE



30/08/2024 09:28:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SPO2 REUSABLE DIGITAL SENSOR FOR ADULT

SUNTECH

Ref. DSU250900

Description

SpO2 articulated finger sensor for adults with cable with a length of 3 meters for Suntech 247 station

Features

| | | | |
|------------------------------------|-----------------|---------------------------|-----------|
| CE Marking | CE0413 | Medical device | Class IIb |
| Bar code | 010869887032265 | HS Code | 90181990 |
| Country | Turkey | | |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 0.052 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

 30/08/2024 09:28:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr