

## 10'3M PATIENT BP HOSE BAYONET TO BAYONET

# SUNTECH


Ref. DSU551456



 30/08/2024 10:12:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## 10'/3M PATIENT BP HOSE BAYONET TO BAYONET

SUNTECH

Ref. DSU551456

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446032562	HS Code	3917.23.0000
Warranty (months)	12		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	2	Packed unit / Height (cm)	18
Packed unit / Length (cm)	19	Packed unit / Weight	0.1
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 30/08/2024 10:12:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)